

## Instructions and User Guide

### Midscapular Thoracentesis Ultrasound Training Model

BPTT2-1005



# Contents

3	<b>Chapter 1: Overview</b>
3	Giving you the confidence only experience can offer™
3	About Blue Phantom™
3	Blue Phantom™ Warranty
4	<b>Chapter 2: Product Cautions</b>
4	Important Safeguards
5	<b>Chapter 3: Product Information</b>
5	Blue Phantom™ Midscapular Thoracentesis Model
5	Included in this Package
5	Additional Items Required for Training
5	Introduction to Your Training Model
6	Quick Facts about Your Training Model
7	<b>Chapter 4: Utilizing Your Training Model</b>
7	Anatomy of Your Training Model
9	Using Your Training Model
9	Using Needles and Catheters
10	Performing Thoracentesis Procedures
10	Maintaining Optimal Fluid Levels
11	How to Determine if Model Requires Refilling
12	Refilling Your Training Model
13	<b>Chapter 5: Caring for Your Training Model</b>
13	Proper Use and Care
13	Swapping Tissue Inserts
14	Cleaning Your Training Model
14	Storing Your Training Model
14	Blue Phantom™ Customer Support

# Overview

# 1

## Giving you the confidence only experience can offer™

Congratulations on the purchase of your Blue Phantom™ ultrasound model(s) for hands-on training. Every product we manufacture at Blue Phantom™ is specifically designed to be the most realistic and ultra-durable ultrasound simulation phantoms available anywhere. Our high standards for quality manufacturing and design guarantee that you receive only the absolute best.

## About Blue Phantom™

Blue Phantom™ brings you the most realistic and durable hands-on ultrasound training models available anywhere. At Blue Phantom™ we know that learning to use ultrasound requires practice. You gain confidence and skill through experience. That is why we offer you the best ultrasound simulation training available.

## Blue Phantom™ Warranty

Blue Phantom™ takes pride in its quality design and manufacturing standards. Our products are warranted to you by Blue Phantom™ for one year from the date of purchase against defects in workmanship and materials. During the warranty period, a defective part or product will be replaced either with a new or reconditioned part or product, depending on the availability at the time.

This warranty covers normal consumer usage and does not cover damage incurred through use not consistent with the product design. Failure that results from alteration, accident, misuse, vandalism, or neglect is not covered under this warranty. This warranty does not extend to any products that have been used in violation of written instructions.

# Product Cautions

# 2

Please read this instruction guide carefully. Do not begin using this model until you fully understand these safeguards and have read the user guide in its entirety.

## Important Safeguards



1. Read Instructions – All safety and operating instructions should be read before the unit is operated.
2. While all parts of this user guide are important, the red flag that you see to the left denotes especially important content. Please familiarize yourself with all of the content prior to using your training model or damage to the model can occur.
3. This user guide includes instructions applicable to a variety of model configurations.
4. Retain Instructions – The safety and operating instructions should be retained for future reference.
5. Heed Warnings – All warnings in the operating instructions should be adhered to.
6. Follow Instructions – All operating and maintenance instructions should be followed.
7. Weight Warning—Use caution as this model is heavy. Use proper lifting techniques to prevent bodily injury. Notify others of the heavy weight warning to prevent others from being injured while operating or moving the model. Product weight is approximately 23lbs (11kg).
8. Care must be taken to place the model in a position in which it will not fall off of the bed or surface, as this may cause injury.
9. Accessories – Do not place this unit on an unstable cart, stand, tripod, bracket, or table. The unit may fall causing serious injury to a child or adult, and serious injury to the unit.
10. **CAUTION:** Please use extreme care when using needles and sharp objects as to not accidentally injure yourself during training.



## Product Information

# 3

### Blue Phantom™ Midscapular Thoracentesis Model

Blue Phantom™ Midscapular Thoracentesis Ultrasound Training Model ([BPTT2-1005](#))

#### Included in this Package

- Blue Phantom™ Midscapular Thoracentesis ultrasound training model ([BPTT2-1005](#))
- Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#))
- User Guide and Utilities CD

#### Additional Items Required for Training

- For optimal performance, please use 18-21 gauge sharp and unbent needles and associated catheter kits (please see *Chapter 4: Utilizing Your Training Model—Using Needles* section on page 9 of this user guide)
- Ultrasound system configured with an appropriate transducer for ultrasound guided thoracentesis procedures
- Ultrasound gel

#### Optional Accessories for Your Training Model

- Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#))
- Soft Storage Case ([BPTT-1011](#))
- Blue Phantom™ Midscapular Thoracentesis Replacement Tissue Insert ([BPTT2-1006](#))

#### Introduction to Your Training Model

This model is intended as a platform for ultrasound guided midscapular thoracentesis hands-on technique training. The model is designed to be extremely realistic, mimicking the feel and imaging characteristics associated with midscapular approach to ultrasound guided thoracentesis procedures. Blue Phantom's self healing design provides you with superb durability. In order to get the most out of your training platform, it is important that you properly care for your model.



*(continued on next page)*

## Quick Facts about Your Training Model

- Excellent platform for training clinicians in the psych-motor skills associated with ultrasound guided thoracentesis procedures
- Ultrasound tissue module contains the chest wall superficial tissue, 6th-9th ribs and intercostal spaces, pleural cavity with lung and atelectatic lung, diaphragm and superior spleen
- Positive fluid offers users immediate feedback when pleural effusion fluid is accurately accessed
- Superb ultrasound imaging characteristics; experience the same quality you expect from imaging patients in a clinical environment
- Made in USA



# Utilizing Your Training Model

# 4

## Blue Phantom™ Midscapular Thoracentesis Model

BPTT2-1005

### Anatomy of Your Training Model

1. Remove your training model from its shipping container and make sure that you have received all of the items listed in *Chapter 3: Product Information—Included in this Package section* on page 5 of this user guide. If you did not receive one of the listed items, or if you received the wrong items, please contact Blue Phantom™ Customer Support immediately:

Telephone: (425)881-8830

Email: [customersupport@bluephantom.com](mailto:customersupport@bluephantom.com)

Web: [www.bluephantom.com](http://www.bluephantom.com)

2. Familiarize yourself with your training model.



Blue Phantom™ ultrasound training models are constructed using our patented simulated ultrasound tissue and mimics imaging characteristics to that of human tissue. Care must be taken to not place the model on rough surfaces as the model can take on the characteristics of that surface. Do not place objects under the model as the tissue is soft and will conform to the shape of the object.

*(continued on next page)*

## Anatomy of Your Training Model (*continued*)

Your Blue Phantom™ Midscapular Thoracentesis model has one tube exiting from the model. The tube acts as a quick fill port for your model and comes pre-filled with Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#)) to represent pleural fluid.

Please see *Image A* below for the location of the quick fill tube.

Quick Fill Tube for  
simulated pleural fluid  
([BRS183-Yellow](#))



*Image A*



Your training model comes pre-filled with Blue Phantom™ Ultrasound Refill Solution. If you suspect that your model is low on fluid, please refer to *Chapter 4: Utilizing Your Training Model—Maintaining Optimal Fluid Levels* section on page 10 of this user guide.



When refilling your model, take care not to lose the blue end cap.

## Using Your Training Model

1. Remove your Blue Phantom™ from its packaging and place on a clean, hard, flat surface. Position your model upright.
2. This model is designed for ultrasound guided technique training.
3. If utilizing ultrasound guidance, place ultrasound gel on the model or on the ultrasound transducer in adequate quantities so that the probe slides effortlessly across the surface of the model. Add more gel as necessary.
4. Adjust the ultrasound system controls per the manufacturer's instructions, increasing and decreasing the depth and gain controls until the desired image is obtained.



If you intend on guiding sharp objects into the phantom, never place the model in a location where you might accidentally puncture yourself.

## Using Needles and Catheters



1. For best performance, we recommend that you utilize standard thoracentesis kits and new, sharp, unbent 18-21 gauge needles and associated catheter kits when accessing the structures in the model.
2. Do not use any needle larger than 18 gauge or permanent damage to your model may occur.
3. Smaller bore needles (>21 gauge) can bend during use and damage your model's simulated tissue.
4. Aggressive repositioning of needles rather than removing and repositioning can cause stubborn or permanent needle tracks due to the needle tip dragging through the simulated tissue.
5. Dull needles may also cause permanent damage to the tissue. It is important to replace needles approximately every ten cannulations.
6. If you strike bone during training, remove the needle and replace it before re-inserting. Striking the bone can cause the needle to become dull.

## Performing Thoracentesis Procedures

The Blue Phantom™ Midscapular Thoracentesis model accommodates full ultrasound guided thoracentesis procedural training, including:

- Create a sterile field
- Introduce needles to the pleural space
- Withdraw simulated pleural fluid - Blue Phantom Yellow Ultrasound Refill Solution ([BRS183-Yellow](#))
- Experience positive fluid flow when pleural effusions are accurately accessed
- Thread small catheters (18-21 gauge needles and associated catheter kits)

Please always utilize small gauge thoracentesis catheter kits and sharp and unbent 18-21 gauge needles. For more information, please refer to *Chapter 4: Utilizing Your Training Model—Using Needles and Catheters* section on page 9 of this user guide. Always heed the following warnings while utilizing your training model:



1. DO NOT use antiseptics such as iodine on your training model. This may cause permanent damage to the simulated tissue.
2. Users can remove fluid from the pleural space when the pleural effusion is accurately cannulated.
3. Please note that any fluid withdrawn from the model will require refilling. Lack of fluid will result in diminished ultrasound image. For instructions on refilling your model, please refer to *Chapter 4: Utilizing Your Training Model—Refilling Your Training Model* on page 12 of this user guide.

## Maintaining Optimal Fluid Levels

The Blue Phantom™ Ultrasound Refill Solution contained within the model is specially formulated to allow for optimal performance. It is important to maintain a good fluid level within your ultrasound training model. Through normal repeated use, users can experience a reduction in the amount of fluid, resulting in small amounts of air within the model. This is directly dependent upon the amount of fluid removed during simulation training.



**WARNING:** Using fluid other than that supplied by Blue Phantom™ will void your warranty. Using other fluid will cause changes in the imaging qualities of the model, reduce the ability to thread catheters, and cause fungal or bacterial growth within the model. DO NOT dilute the Ultrasound Refill Solution contained within your model.

Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#)) is utilized within the model to simulate pleural fluid.



The simplest way to maintain optimal fluid level in the vessels is to have users inject the accessed fluid back into the model after accessing the pleural effusion.

*(continued on next page)*

## Maintaining Optimal Fluid Levels (*continued*)



**CAUTION:** Use refill solution only as directed. Not intended for human consumption. If accidental consumption occurs, drink a glass of water and consult a physician. May irritate eyes; flush well with water. May contain pigments that can stain clothing; wash immediately with cold soapy water. Keep out of reach of children.


Please see below to determine if your model requires refilling and for complete refilling directions.

### How to Determine if Your Model Requires Refilling


Choose one of the following methods to determine if the pleural cavity requires refilling:

1. **Ultrasound Method**

An optimally filled pleural cavity will be identified by the presence of a black echo-free lumen (refer to *Image A*). A low fluid environment is identified by the inability to visualize the vessels during normal imaging situations (refer to *Image B*). This is due to the presence of air within the vessels, which will reflect all of the sound energy.



*Image A*



*Image B*

2. **Non-Imaging Method**

The presence of air in the Quick Fill Tube and lack of pleural fluid flow are indicators that your model requires refilling. Hold the Quick Fill Tube upright towards the ceiling in the fully extended position for at least 60 seconds. If an air bubble migrates to the terminal end of the tube, this indicates that the model requires refilling.

If you have determined that your model requires refilling, please refer to *Chapter 4: Utilizing Your Training Model—Refilling Your Training Model* section on page 12 of this user guide.

## Refilling Your Training Model

Refilling your Blue Phantom™ Midscapular Thoracentesis ultrasound training model is a simple process that will take approximately 10 minutes.



If you are uncomfortable refilling your training model for any reason, contact Blue Phantom™ and we can refill the phantom for a nominal fee.

### Items Required

- Empty syringe
- IV Bag (*optional*)
- Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#))

### Directions for Refilling Your Training Model

1. Position the training model in the upright position.
2. Fill the empty syringe with Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#)).
3. Remove the blue end cap from the Quick Fill Tube and hold the Quick Fill Tube upright towards the ceiling in the fully extended position. Hold the tube upright for at least 60 seconds to allow any accumulated air to migrate to the terminal end of the tube.
4. Connect the filled syringe's luer lock female connector to the male connector on the Quick Fill Tube.
5. SLOWLY infuse the refill solution into the Quick Fill Tube in 10ml increments until it is full. Be sure to purge air after each successive 10ml by pulling back on the plunger. Filling the Quick Fill Tube SLOWLY prevents air bubbles from being introduced into the model.
6. High volume users will benefit from connecting either a syringe or an IV bag containing Blue Phantom™ Ultrasound Refill Solution to the Quick Fill Tube. As users withdraw fluid from the model, the Quick Fill Tube is automatically refilled. To do so, please refer to steps 8-10 below. Otherwise, replace the blue cap back on the luer lock connector.
7. To easily maintain optimal fluid levels, fill a syringe or an empty IV bag with Blue Phantom™ Yellow Ultrasound Refill Solution ([BRS183-Yellow](#)).
8. Connect the full syringe or IV bag to the male luer lock connector on the Quick Fill Tube.
9. If utilizing an IV bag, be sure to hang the IV bag 1ft above the model.



When refilling the Quick Fill Tube or utilizing a full syringe or IV bag, take care not to lose the blue end cap.

# Caring for Your Training Model

# 5

## Proper Use and Care

Proper care of your training model will result in tremendous utility. Please heed all instructions contained in this user guide when using your model.

## Swapping Tissue Inserts

Swapping tissue inserts in and out of your Blue Phantom™ Thoracentesis ultrasound training model is a simple process that will take approximately 10 minutes.

### Items Required

- Empty syringe
- Soapy water mixture: 1ml liquid soap, 4ml tap water
- Paper towels

### Directions

1. Position your Midscapular Thoracentesis training model on a clear workspace in the upright position.
2. Fill the empty syringe with the soapy water mixture.
3. Inject the soapy water solution into the holes where the Quick Fill Tube exits from the base platform of the training model. This will lubricate the tube for easy removal.



DO NOT inject the soapy water mixture into the Quick Fill Tube.

4. Holding onto the top of the tissue insert, gently lift the tissue insert module superiorly so that it pulls both the tissue insert and the Quick Fill Tube out of the model.
5. Completely remove the tissue insert module and the Quick Fill from the base platform.
6. Now you are ready to insert the new tissue insert. Please refer to *Image A* on page 8 of this user guide to determine the proper orientation of the tissue insert.
7. Guide the Quick Fill Tube into the hole located in the recess in base platform. Thread the tubes until the tissue insert reaches the recess in the base platform. Gently position the tissue insert into the recess until it is positioned flush with the surface of the thoracentesis base platform.
8. Make any necessary adjustment to the tissue insert and Quick Fill Tube so that the insert and tube are properly positioned.



For more information, please refer to our instructional YouTube movie on [How To Replace Ultrasound Tissue Inserts](#).

## Cleaning Your Training Model

After each use, your training model can be easily cleaned using mild soapy water. For best results, mix one part liquid soap with one part tap water. Gently rinse the model with the soapy water to remove any accumulated debris.

Use a clean, soft, lint-free cloth to dry after cleaning. Dry the model using a dabbing motion, rather than wiping or rubbing the model.



Wiping or rubbing the surface aggressively can result in scuffing the simulated tissue.

## Storing Your Training Model

The model can be stored at room temperature either in the open or in the Blue Phantom™ soft storage case (optional). Always store your Midscapular Thoracentesis training model in the upright position.



Do not store the model in contact with other objects. This can cause the simulated tissue to become deformed.

## Blue Phantom™ Customer Support

Blue Phantom™ is committed to providing you with superb products and uncompromising customer support. Should you require assistance feel free to contact us directly at:

Telephone: (425)881-8830

Email: [customersupport@bluephantom.com](mailto:customersupport@bluephantom.com)

Web: [www.bluephantom.com](http://www.bluephantom.com)



Find us on  
**Facebook**

