CAE CathLabVR™ User Guide Administrator's Guide



905K400552- v2.0 April 2021

PROPRIETARY NOTICE: This document, including the information contained herein, is confidential and/or proprietary to CAE Inc., and shall not be reproduced or disclosed in whole or in part, or used for any purpose whatsoever without the prior written authorization of CAE Inc.



Your worldwide training partner of choice

Proprietary Notice

This document, including the information contained herein, is confidential and/or proprietary to CAE Inc., and shall not be reproduced or disclosed in whole or in part, or used for any purpose whatsoever without the prior written authorization of CAE Inc.

Trademark Information

Trademarks and/or registered trademarks of CAE Inc. and/or its affiliates include but are not limited to Apollo, Ares, AresAR, Athena, BabySIM, Caesar, CathLabVR, EndoVR, HPS Human Patient Simulator, iStan, Juno, LapVR, LearningSpace, Lucina, LucinaAR, Luna, Maestro, Muse, NeuroVR, Orion, PediaSIM, Vimedix, VimedixAR and Vivo. All other brands and product names are trademarks or registered trademarks of their respective owners. All logos, tradenames and trademarks referred to and used herein remain the property of their respective owners and may not be used, changed, copied, altered, or quoted without the written consent of the respective owner. All rights reserved.

CAE ©2021



CONTENTS

Introduction
Contained in this Guide2
CathLabVR Software2
Cautions and Warnings
General
Operations
Ergonomics
Using the Administrator Tools5
Starting the Software5
Interface Overview
Hiding the Side Panel
Users and Groups
Creating Users
Editing User Information
Deactivating Users
Reactivating Users
Modifying Groups
Emailing Users and Groups24
Curriculum
Creating a Course
Editing a Course
Creating a Task
Renaming a Task
Assigning a User to a Task42
Assigning a Group to a Task43
Assigning a Course to a Task44
Deleting a Task
Faculty Content
Configuration
Modifying the Network Settings
Backing Up the Database50
Restoring the Database
Adjusting the Volume Controls

Administrator's Guide

Print
Print Preview
Printing Data to a PDF File60
Printing Data to a Network Printer62
Import/Export
Importing Users and Groups64
Exporting Users and Groups67
Exporting Results
Exiting the Software
Appendix A - FAQ
Dual Monitor Configuration



INTRODUCTION

The CAE CathLabVR[™] simulator is one of a series of didactic tools that can be used to improve training for endovascular diagnostic and intervention procedures. Using modified endovascular tools, trainers can navigate through virtual anatomy based on actual patient data. The simulator is not intended to replace any courses or hands-on sessions with live subjects, but should be implemented as part of a blended learning solution.

IMPORTANT: The simulator is not intended to diagnose the condition of a live human being or identify a life-threatening situation.

The simulators include:

- **CathLabVR hardware interface device** This device is instrumented to track the manipulation of coaxial endovascular tools and provides computer-controlled feedback based on the interaction with simulated patient anatomy.
- Laptop computer for running the simulation software with a Bluetooth[®] mouse.
- **Modified endovascular tools**, such as catheters, leads and guide wires for performing different simulated endovascular procedures.



The CathLabVR Simulator

Contained in this Guide

This Administrator Guide has been designed for quick access to information on how to use the Administrative Tools available in the CathLabVR software.

- The **Cautions and Warnings** section is for the safety of users as well as for the protection of the simulator. Please be sure to read this prior to using the simulator. Follow the guidance during use.
- The **Using the Administrator Tools** section provides instructions on the use of the various administrative software features such as, creating and maintaining users and groups, and creating curriculum.
- Appendix A: **FAQ** provides information regarding common issues involving simulation activities. t

Note: Refer to the *CAECathLabVR User Guide* for your specific product regarding information on running simulations and equipment related topics as follows:

- Specifications
- Equipment Overview
- Setup
- Care and Maintenance

CathLabVR Software

The following identifies the software version associated with the release of this guide.

Document	Document Version	Software
CathLabVR Administrator's Guide	905K400552 v2.0	UniSim



CAUTIONS AND WARNINGS

Please read and understand these cautions and warnings before you begin using the simulator.

General

- When unpacking the simulator for the first time, use box cutters carefully to protect both the packaging and the product.
- Do not use equipment in a manner that is not specified by the manufacturer.
- Do not operate the equipment with any covers removed.
- Keep all electronic equipment free from high humidity or liquids.
- Do not leave the hardware powered on when the simulator is not in use.

Operations

- Do not handle the tracking units while in use and never touch the gold metallic parts.
- Always ensure laptop is plugged into a power source during simulator operation. Batterypower is insufficient to run the laptop during simulation.
- Do not use or connect any USB device or cable to the simulator or laptop other than the USB cable provided.
- When installing equipment near an electric socket, ensure cables do not represent a hazard.
- Avoid weighing down the equipment with any tools or accessories not relevant to the immediate exercise.
- Modified endovascular tools included with the simulator can be very sharp. Pay attention while handling them.
- Refrain from bending, kinking, or stepping on the endovascular tools.
- Do not insert an end of an endovascular guidewire into a power outlet.
- During tool manipulation the endovascular guidewire extends out of the rear end module. Ensure that the area around the rear end module is clear of obstacles or people.
- Ensure you have a table with the minimum recommended width and length to support the Tracking System and Laptop.

Ergonomics

- Improper or prolonged keyboard use may result in injury.
- Viewing a monitor screen for extended periods of time may result in eye strain.

THIS PAGE INTENTIONALLY LEFT BLANK



USING THE ADMINISTRATOR TOOLS

Administrators managing the CathLabVR simulator are responsible for setting up user accounts, privileges and curriculum. Additionally, the administrator has the ability to change network settings, set up emails, generate and print reports, and import and export system data.

Note: For optimal use, no other software programs should be open while the simulator software is running.

Starting the Software

1. Once the CathLabVR laptop is powered on, the CAE Start screen displays.



CAE Start Screen

2. Double-click the **CathLabVR** shortcut icon.

3. Accept the End User License Agreement. The login screen appears.


	CAE CathLab Reries
The Username field——	Username
The Password field	Password
The Login button ——	Login
	Copyright © CAE Healthcare Inc., 2019. All Rights Reserved. 1.0.13.0

The Login Screen

To log in to the software:

- 1. Enter the assigned username in the **Username** field. The default username for a new administrator is *admin*.
- 2. Enter the assigned password in the **Password** field. The default password for a new administrator is *admin*.
- 3. Click Login.

The administrator Home screen appears, indicating the administrator is logged in.




The Administrator Home Screen



From the administrator Home screen, administrators can navigate to the screens required for editing user and group accounts, assigning curriculum, adjusting configurations and printing, importing and exporting user data.

Interface Overview

When the Administrator first logs in to the software, the Home screen is displayed and the **Users & Groups** icon is selected by default.



The Administrator Home Screen

The CathLabVR administrator interface contains a dashboard of icons located near the top of the Home screen that represents functions available for administrators to perform various tasks.



The Icon Dashboard

When an icon is selected, the associated content appears in the main panel of the interface. The side panel, on the left of the screen, displays available selections which change based on the icon chosen. For example, when the **Users & Groups** icon is chosen, the side panel is referred to as the *Users panel* and has selections associated with managing users and groups. When the **Configuration** icon is chosen, the side panel is referred to as the *Configuration* icon is chosen, the side panel is referred to as the *Configuration* panel and has selections associated with configuration panel and has selections associated with configuration panel and has selections associated with configuration panel and has selections associated with configuration.

Users & Groups

From the Users & Groups screen, administrators can create user accounts, assign privileges, courses and tasks to individual users and assign users, courses and tasks to groups. Administrators can also deactivate user accounts and set up and send out automated emails containing user account information.

Curriculum

From the Curriculum screen, administrators can assign courses and tasks to individual users and assign courses, tasks and users to groups. Administrators can also upload faculty content to appear in the didactic content of a simulation, create new tasks and courses, and modify user parameters.

Configuration

From the Configuration screen, administrators can restore or backup the database, adjust network settings, and adjust language and volume preferences.

Print

From the Print screen, administrators can print user account information, usage summaries and reports for individual users and groups.

Import/Export

From the Import/Export screen, administrators can import users and export user lists and results.



Hiding the Side Panel

Administrators can hide the side panel to view content in full screen mode.



The Home Screen

To hide the side panel:

1. From the side panel, select the desired item.

The content appears in the main screen.

2. Click the Hide Side Panel arrow.

The content expands to the width of the full main screen.

Clicking on the **Hide Side Panel** arrow again will redisplay the side panel of options.

Users and Groups

The Users and Groups tabs gives access to create, modify, and remove user and group profiles.

Administrators are responsible for creating and managing users and groups in the CathLabVR system. Only users with administrative privileges can create new users and groups.

Creating Users

To create a user:


1. From the icon dashboard, click the **Users & Groups** icon.



The Users & Groups Icon

The Users & Groups screen appears.


- 2. From the Users and Groups screen, click the **Users** tab.
- 3. Click the **Add User** button.



The Users & Groups Screen



The Personal Data screen appears.



The Personal Data Screen


- 4. Enter the information for the specific user in the different fields.
- 5. Click Done.

Assigning Courses to a User Account

To assign a course to a user from the Personal Data screen:

1. Select the user from the **Users** panel on the Users and Groups screen.

The Personal Data screen appears.



The Personal Data Screen

- 2. Click on the Courses & Tasks tab in the Courses & Tasks panel.
- 3. Click on the **Courses** tab.
- 4. Select the desired course in the **Courses & Tasks** panel.
- 5. Hold down the left mouse button and drag the course from the **Courses & Tasks** panel to the **Assigned Courses & Tasks** panel.

Once the button is released, the course appears in the assigned curriculum for the user.

Assigning Tasks to a User Account

To assign a task to a user from the Personal Data screen:

1. Select the user from the **Users** panel on the Users and Groups screen.

1 -**Q**-₽ -The Users & Groups icon The Tasks view The Courses view The Courses & The Assigned-Tasks panel **Courses &** Tasks panel The Courses & Tasks tab The Users The **Groups** tab tab Г 11 The Personal Data Screen

The Personal Data screen appears.

- 2. Click on the **Tasks** tab.
- 3. Select the desired task in the **Courses & Tasks** panel.
- 4. Hold down the left mouse button and drag the task from the **Courses & Tasks** panel to the **Assigned Courses & Tasks** panel.

Once the button is released, the task appears in the assigned curriculum for the user.



Editing User Information


To edit a current user's information:

1. From the icon dashboard, click the **Users & Groups** icon.



The Users & Groups Icon

The Users & Groups screen appears.



The Users & Groups Screen

- 2. From the Users and Groups screen, click the **Users** tab.
- 3. Select the desired user from the **Users** panel.

The Personal Data screen appears.

	CAE Healthcare - CathLabVR™ ▼		1/7/2014 1:51:42 PM		⑦ Default Administrator • en •	
The Users &—	Users & Groups Curriculum	Configuration Print Imp	C			The Edit
Groups icon	USRS O	JohnDoe Personal Data User ID: JohnDoe First Name: John	Password: *** Last Name: Doe	Date of Birth: 2/12/1985 Experience: Medical Student	Edit Print Deactivate COURSES & TASKS O COURSES COURSES Sample Course 1 -	button
The Users —		Email: johndoe@cae.com Administrator: No	Phone: Data Recorded: Yes	Dominant Hand: Right	Sample Task Sample Course 2 Sample Task Sample Course 3 Sample Task	
panel		Assigned Course & Taiks Sample Course 3	Ansigned Groups Part Year Residents	Group Course & Takis First Year Residents Sample Course 3 Task - Case 1	 Sample Taik-Care J Sample Taik-Care A Sample Taik-Care S Sample Taik-Care B Sample Taik-Care B Sample Taik Sample Taik Sample Taik 	
The Users tab	CAE Institutes					


The Personal Data Screen

4. Click on the **Edit** button.

The Personal Data screen is now accessible for editing.

- 5. Change the desired information for the specific user in the different fields.
- 6. Click Done.

Note: If the user navigates away from the Personal Data screen before clicking **Done** to save the information, the **Unsaved Personal Data** message will appear.



The Unsaved Personal Data Message

To save the changes, click **Save**.

To discard the changes, click **Discard**.




Deactivating Users

If a learner is no longer attending training or their account is no longer active, the administrator can deactivate the user.

To deactivate a user:

- 1. From the Users and Groups screen, click the **Users** tab.
- 2. Select the desired user.

The user's information appears in the personal data screen.



The Inactive Users Personal Data Screen

- 3. Click on the **Deactivate** button.
- 4. At the confirmation message, click on the **Deactivate User** button.

The user account is deactivated and the user's name will no longer appear in the User panel. Although, deactivated users may be reactivated; refer to the *Reactivating Users* section of this guide.

Deleting Users

A user may be deleted by an administrator, only after deactivating.

To delete a user:

- 1. If the user is currently active, follow the previous steps to deactivate the user.
- 2. From the Users and Groups screen, click the **Users** tab.
- 3. Click the **Display Inactive Users** button.

The deactivated user names appear in the **Users** panel in gray, italic font.

Note: All other buttons in the Users panel must be toggled off for the delete option to be available. The only active button should be the Display Inactive Users. Should other buttons appear blue, click to toggle off.

- 4. Select the desired deactivated user.
- 5. Click on the **Delete** button.
- 6. At the confirmation window, select the **Delete User** button.

Note: Once the user is deleted, the action cannot be undone.



Reactivating Users

Upon creating and saving a user, the administrator activates the user. If a user has been deactivated by an administrator, they can also be reactivated.

To activate a currently inactive user:

- 1. From the Users and Groups screen, click the **Users** tab.
- 2. Click the **Display Inactive Users** button.




The Users and Groups Screen

The deactivated user names appear in the **Users** panel in gray, italic font.

3. Select the desired deactivated user.

The user's information appears in the personal data screen.



The Inactive Users Personal Data Screen

4. Click Activate.

Creating Groups

To create a group:

1. From the icon dashboard, click the **Users & Groups** icon.




The Users & Groups Icon

The Users & Groups screen appears.




2. From the Users and Groups screen, click the **Groups** tab.



The Users & Groups Screen - Groups Tab

3. Click the **Add Group** button.

The New Group appears in the **Groups** panel.



The Edit Group Screen


- 4. Enter the desired group name in the **Group Name** field.
- 5. Click **Done** to save the group.

Assigning a User to a Group

To assign a user to a group from the Personal Data screen:

1. Select the user from the **Users** panel on the Users and Groups screen.

The Personal Data screen appears.



The Personal Data Screen

- 2. Click on the **Groups** tab in the **Groups** panel.
- 3. Select the desired group in the **Groups** panel.
- 4. Hold down the left mouse button and drag the group from the **Groups** panel to the **Assigned Groups** panel.

Once the button is released, the group appears in the assigned groups for the user.



Modifying Groups

Administrators can modify the group name, assigned curriculum and assigned users.

To modify a group name:

- 1. From the Users and Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel. The group information screen appears.
- 3. Click Rename.
- 4. Enter the new group name in the **Group Name** field.
- 5. Click Done.

Deleting a Group

To delete a group, click the **Delete** button on the group information screen.

Adding a Course to the Group Curriculum

To add a course to the group curriculum:

- 1. From the Users & Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel.
- 3. From the Group Information screen, click the **Courses & Tasks** tab in the **Courses & Tasks** panel.
- 4. Click on the **Courses** view and select the desired course in the **Courses & Tasks** panel.
- 5. Hold down the left mouse button and drag the course from the **Courses & Tasks** panel to the **Assigned Courses & Tasks** panel.

Once the button is released, the course appears in the assigned curriculum for the group.

Deleting a Course from the Group Curriculum

To delete a course from the group curriculum:

- 1. From the Users & Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel.
- 3. From the **Assigned Courses & Tasks** panel, select the desired course.
- 4. Click the trash can icon button in the upper right corner of the **Assigned Courses & Tasks** panel.

Adding a Task to the Group Curriculum

To add a task to the group curriculum:

- 1. From the Users & Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel.
- 3. From the group information screen, click on the **Courses & Tasks** tab in the **Courses & Task**s panel.
- 4. Click on the Tasks view and select the desired task in the Courses & Tasks panel.
- 5. Hold down the left mouse button and drag the task from the **Courses & Tasks** panel to the **Assigned Courses & Tasks** panel.

Once the button is released, the task appears in the assigned curriculum for the group.

Deleting a Task from the Group Curriculum

To delete a task from the group curriculum:

- 1. From the Users & Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel.
- 3. From the **Assigned Courses & Tasks** panel, select the desired task.
- 4. Click the trash can icon button in the upper right corner of the **Assigned Courses & Tasks** panel.

Adding a User to a Group

To add an individual user to the group:

- 1. From the Users & Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel.
- 3. From the group information screen, click on the **Users** tab in the **Users** panel.
- 4. From the **Users** panel, select the desired user.
- 5. Hold down the left mouse button and drag the user from the **Users** panel to the **Assigned Users** panel.

Once the button is released, the user appears as an assigned user to the group.



Deleting a User from the Group

To delete a user from the group:

- 1. From the Users & Groups screen, click the **Groups** tab.
- 2. Select a group from the **Groups** panel.
- 3. From the **Assigned Users** panel, select the desired user.
- 4. Click the trash can icon button in the upper right corner of the **Assigned Users** panel.

Emailing Users and Groups

Administrators have the ability to coordinate automated email correspondence with individual users and groups. Administrators can customize email settings to include content such as user profiles, reports and usage summaries, provide a return email address and set up a recurring auto update.

Once the email settings are customized, the administrator can send the email immediately or save the email to send at a later time.

Note: To send emails from the simulator, a network connection is required and the administrator must configure the SMTP network information in the simulator. For more information about SMTP configuration, refer to *Adjusting the Network Settings*.

Emailing Users

Administrators can send emails containing user profiles, usage summaries and reports to individual users.

To send an email to an individual user:

1. From the icon dashboard, click the **Users & Groups** icon.



The Users & Groups Icon

The Users & Groups screen appears.



2. From the Users and Groups screen, on the **Users** tab, click the **Emails** button.



The Users & Groups Screen

The Email Options screen appears.

	CAE Healthcare - CathLabVR™ ▼		1/8/2014 1:32:09 PM		⑦ Default Administrator • en •	
	Users & Groups Curriculum	Configuration	Import/Suport			
	USERS E 📼 🗅 🕼 🕊 admin	Choose Send Email Options All Users			Back to User Data	
	guest JohnDoe	Send to: Users Administrator(s) Other:		Content: User Profile Summary Usage Reports		
		Return Address:				
		E Set up a recurring auto updat				
The						The Send Now button
Users tab	AE matricer				(Send Now) (Save	The Save —button

The Email Options Screen

Note: If a user is selected from the **Users** panel, the Email Options screen displays a tab specifically for the user. Otherwise, the **All Users** tab displays.

- 3. Select the recipients for the email.
- 4. Select the content for the email. Content options include **User Profile**, **Summary Usage** and **Reports**.
- 5. Enter a return email address in the **Return Address** text box, if desired.

Using the Administrator Tools

- 6. Select the Set a recurring auto update checkbox, if desired.
- 7. Click **Send Now** to send the email immediately.

Note: Email settings can be customized and saved to send at a later time. To save an email to send later, click **Save**.

Emailing Groups

Administrators can send emails containing user profiles, usage summaries, and reports to groups.

To send an email to a group:


1. From the icon dashboard, click the **Users & Groups** icon.



The Users & Groups Icon

The Users & Groups screen appears.

2. From the Users & Groups screen, on the **Groups** tab, click the **Emails** button.



The Users & Groups Screen - Groups Tab



The Email Options screen appears.

CAE Healthcare - CathLabVR™ ▼		1/8/2014 1:33:45 PM	② Default Administrator • en •	
	🌣 📥 🛟			
	Configuration Print Import/Export			
GROUPS C	Choose Send Email Options		Back to Group Data	
First Year Residents	All Users			
	Send to:	Content:		
	🖾 Users	🗉 User Profile		
	C Administrator(s)	Summary Usage		
	Cther:	III Reports		
	Return Address:			
	E Set up a recurring auto update			
				The Send
				Now button
				NOW DULLON
~				The Save
CAE methors				
1 1			Send Now Save	—button

The Email Options Screen

- 3. Select the group name tab or **All Users** tab.
- 4. Select the recipients for the email.
- 5. Select the content for the email. Content options include **User Profile**, **Summary Usage** and **Reports.**
- 6. Enter a return email address in the **Return Address** text box, if desired.
- 7. Select the **Set a recurring auto update** checkbox, if desired.
- 8. Click Send Now to send emails immediately

OR

Click **Save** to retain the chosen options and send emails at a later time. Click **OK** at the **Data Saved Successfully** screen.

9. Click **Back to User Data** to return to the User panel.

To edit saved email options or to send emails using previously saved options:

1. From the Users & Groups screen, on the **Groups** tab, click the **Emails** button.

The Email Options screen will appear.

2. Edit options and click Save

OR

Click **Send Now** to send emails.

Curriculum

The Curriculum screen enables administrators to create, modify, and assign curriculum courses and tasks.

Creating a Course

To create a course:

1. From the icon dashboard, click the **Curriculum** icon.



The Curriculum Icon


The Curriculum screen appears.



The Curriculum Screen



2. From the Curriculum screen, click the **Courses** tab.



The Course Information Screen

3. Click the **Add Course** button.

The Course Information screen appears.

- 4. Enter the course name in the **Course Name** field.
- 5. If desired, select the **Mandatory Order** checkbox to require assigned users to complete the tasks in the order specified for the course.
- 6. Click Done.

The course is created.

Editing a Course

Once a course has been created, the administrator can modify any of its features, including assigning or removing tasks, users or groups as well as changing the course name.

Renaming a Course


To rename a course:

1. From the icon dashboard, click the **Curriculum** icon.



The Curriculum Icon

The Curriculum screen appears.




The Curriculum Screen

- 2. From the Curriculum screen, click the **Courses** tab.
- 3. Select a course from the **Courses** panel.



The Course Information screen appears.



The Course Information Screen

- 4. Click on the **Rename** button.
- 5. Enter the new course name in the **Course Name** field.
- 6. If desired, select the **Mandatory Order** checkbox to require assigned users to complete the tasks in the order specified for the course.
- 7. Click Done.

The course name is changed and will appear with the new name in previous associations (e.g. user course assignments).

Deleting a Course

To delete a course, click the **Delete** button on the Course Information screen.

The course is deleted.


Assigning a Task to a Course

To assign a task to a course:

- 1. From the Curriculum screen, click the **Courses** tab.
- 2. Select a course from the **Courses** panel.

The Course Information screen appears.

3. From the Course Information screen, click on the Tasks tab in the Tasks panel.



The Course Information Screen

- 4. Select the desired course in the **Tasks** panel.
- 5. Hold down the left mouse button and drag the course from the **Tasks** panel to the **Assigned Tasks** panel.

Once the button is released, the task appears in the assigned curriculum for the course.

Deleting a Task from a Course

To delete a task from a course:

- 1. From the **Assigned Tasks** panel, select the desired task.
- 2. Click the trash can icon button in the upper-right corner of the **Assigned Tasks** panel.

The task is deleted.

Assigning a User to a Course


To assign a user to a course:

- 1. From the Curriculum screen, click the **Courses** tab.
- 2. Select a course from the **Courses** panel.

The Course Information screen appears.



3. From the Course Information screen, click on the **Users** tab in the Users panel.



The Course Information Screen

- 4. Select the desired user in the **Users** panel.
- 5. Hold down the left mouse button and drag the task from the Users panel to the **Assigned Users** panel.

Once the button is released, the task appears in the assigned curriculum for the course.

Deleting a User from a Course

To delete a user from a course curriculum:

- 1. From the **Assigned Users** panel, select the desired user.
- 2. Click the trash can icon button in the upper-right corner of the **Assigned Users** panel.

The user is deleted.


Assigning a Group to a Course

To assign a group to a course:

- 1. From the Curriculum screen, click the **Courses** tab.
- 2. Select a course from the **Courses** panel.

The Course Information screen appears.

3. From the Course Information screen, click on the **Groups** tab in the **Groups** panel.



The Course Information Screen

- 4. From the **Groups** panel, select the desired group.
- 5. Hold down the left mouse button and drag the group from the **Groups** panel to the **Assigned Groups** panel.

Once the button is released, the group appears as an assigned group to the course.

Deleting a Group from a Course

To delete a group from the course curriculum:

- 1. From the **Assigned Groups** panel, select the desired group.
- 2. Click the trash can icon button in the upper-right corner of the **Assigned Groups** panel.

The group is deleted.



Creating a Task

The simulator's module packages come with pre-programmed tasks for each type of procedure. Additionally, administrators can create custom tasks from the **Tasks** tab or the **Task Templates** tab on the Curriculum screen. These tasks are known as administrator-created tasks.


To create a task from the **Tasks** tab:

1. From the icon dashboard, click the **Curriculum** icon.



The Curriculum Icon

The Curriculum screen appears.




The Curriculum Screen

2. Click the **Tasks** tab.


Using the Administrator Tools

3. Click the **Add Task** button.



The Curriculum Screen

The Add Task box appears.



The Add Task Screen

- 4. Select a task to add.
- 5. Click Select.



The new task appears in the **Tasks** panel and the new task information appears in the main panel on the Curriculum screen.

6. Click on the **Done** button.

Note: Duplicate tasks names are not permitted.


To create a task from the **Task Templates** tab:

1. From the icon dashboard, click the **Curriculum** icon.



The Curriculum Icon

The Curriculum screen appears.



The Curriculum Screen

2. From the Curriculum screen, click the **Task Templates** tab.

The Task Template screen appears.

	CAE Healthcare - CathLabVR™ ▼		1/7/2014 9:31:30 AM		⑦ Default Administrator • en •
	Users & Groups Curriculum	Configuration	2000rt/Sport		
	TASK TEMPLATES	Sample Course 1			Create
A selected	* Sample Course 1	User Parameters		Target Performance	
task_	Sample Task Case 1			*	
cubit_	Case 1 Case 2	2			
	Case 3	>			
The Tasks —	Case 4	2			
Templates	 Sample Task Sample Course 2 				
remplates	Sample Task				
panel	 Sample Course 3 Sample Task 				
•	Sample Task				
	Sample Course 4				
	Sample Task Sample Task				
	Sample task				
	CAE Human				
	2 9 0		C	ancel	

The Task Template Screen

3. From the Task Template screen, select a task in the Task Templates panel.

The task parameters for the selected task appear on the Task Template screen.

4. Click **Create**.

The Create Task box appears.

Create Task		
Task Name	Sample Task	
Case Name	Case 1: 57 year-old MALE	
		—The Create
	Create Task Cancel	Task button

The Create Task Box

5. Enter the desired task and case name information for the new task.

Note: New tasks that are not renamed, will appear with the original task template name and "Copy ###" in the **Task Templates** panel. This task name can be edited at a later time in the **Tasks** tab after it is created.


6. Click Create Task.

The Task Template screen appears.

Note: The new task does not appear in the Task Templates panel.



7. From the Task Template screen, click the **Tasks** tab.



The Tasks Screen


The new task appears in the **Tasks** panel on the Tasks screen.

Renaming a Task

Tasks created by the administrator or tasks predefined by the system may be assigned a new name.

To rename a task:

1. From the **Tasks** panel, select the desired task.




The Tasks Screen

The task information appears in the main panel of the Edited Tasks screen.

2. Click on the **Rename** button.



The task information fields become editable.



The Tasks Screen

- 3. Enter the desired named for the new task in the **Task Name** text box.
- 4. Enter the desired case name in the **Case Name** text box.
- 5. Click on the **Done** button.


The task name is updated in the **Tasks** panel.

Assigning a User to a Task

To assign a user to a task:

1. From the **Tasks** panel, select the desired task.

The task information appears on the Tasks screen.



The Tasks Screen

2. From the Tasks screen, click the **Users** tab.

A list of available users appears.

- 3. Select the assigned user in the **Users** panel.
- 4. Drag the user by holding down the left mouse button and navigating the cursor to the **Assigned Users** panel.

The task is assigned to the user.






Assigning a Group to a Task

To assign a group to a task:

1. From the **Tasks** panel, select the desired task.

The task information appears on the Tasks screen.



2. From the Tasks screen, click the **Groups** tab.

A list of available groups appears.

- 3. Select the desired group in the **Groups** panel.
- 4. Drag the group by holding down the left mouse button and navigating the cursor to the **Assigned Groups** panel.

The task is assigned to the group.


Assigning a Course to a Task

Note: Tasks are assigned to courses automatically based on the task template used to create them. Administrators can modify these course assignments for administrator-created tasks.

To assign a course to a task:

1. From the **Tasks** panel, select the desired task.

The task information appears on the Tasks screen.



The Tasks Screen

- 2. Click the **Courses** tab.
- 3. Select the desired course in the **Courses** panel.
- 4. Drag the course by holding down the left mouse button and navigating the cursor to the **Assigned Courses** panel.

The task is assigned to the course.



Deleting a Task

Administrators can delete tasks that are not pre-programmed in the modules (i.e., administratorcreated tasks).

To delete a task:

1. From the **Tasks** panel, select the desired task.




The Tasks Screen

The task information appears in the main panel of the Edited Tasks screen.

2. Click on the **Delete** button.

The Delete Task confirmation message appears.



The Delete Task Message

3. Click on the **Delete Task** button to delete the task.

The task is deleted and no longer appears in the **Tasks** panel.

Faculty Content

Administrators can upload content that will appear on a separate tab labeled **Faculty Content** within the didactic content options.

Adding Faculty Content

To upload documents for the **Faculty Content** tab:


1. From the icon dashboard, click the **Curriculum** icon.



The Curriculum Icon

The Curriculum screen appears.

1. From the Curriculum screen, on the **Tasks** tab, select the desired task.



The Curriculum Screen

2. Click the **Faculty Content** button.



The Select Faculty Content File dialog box appears.

	Select Faculty Content File	
	Faculty Content File Path Browse	—The Browse button
The Update — button	Update Remove Cancel	

The Select Faculty Content File Dialog Box

3. Click **Browse**.

The **Select File** dialog box appears.

Note: If the file is located on an external device, use the USB port located on the back of the simulator to connect the external device to the simulator.

- 4. Select the .htm or .html file to be uploaded.
- 5. From the **Select File** dialog box, click **Open**.

The file name appears in the **Select Faculty Content File Path** field of the **Select Faculty Content File** dialog box.

6. From the **Select Faculty Content File** dialog box, click on the **Update** button.

When users access the task, the **Faculty Content** tab appears and the content from the chosen file populates the screen.

Removing Faculty Content

To remove Faculty Content associated with a task:

- 1. Repeat steps 1-5 as detailed in Adding Faculty Content.
- 2. From the **Select Faculty Content File** dialog box, click on the **Remove** button.

When users access the task, the content will no longer be associated with the task.

Configuration

The Configuration tab gives access to settings controls for the simulator. Through this tab, the administrator can modify the network settings, manage the database, choose language preferences, and adjust volume settings.

Modifying the Network Settings

The Network Configuration screen is used to set the email server address. The server name must be for an outgoing (SMTP) mail server. This server should provide access to all users and administrators who are to receive reports through the email function on the Users and Groups screen.


To modify the network settings:

1. From the icon dashboard, click the **Configuration** icon.



The Configuration Icon


The Configuration screen opens.



The Configuration Screen



2. From the Configuration screen, select **Network Configuration** in the **Configuration** panel.



The Network Configuration Screen

- 3. Enter the SMTP server address in the **Server address** field.
- 4. Select a database option.
- 5. Click **Restart System** to save the information and update the settings.

Backing Up the Database

Administrators can back up the simulator database using the Database Backup feature located on the Configuration screen.

To back up the simulator database:


1. From the icon dashboard, click the **Configuration** icon.



The Configuration Icon

2. From the Configuration screen, select the **Database Backup/Restore** option from the **Configuration** panel.

The Database Backup/Restore screen appears.



The Database Backup/Restore Screen

3. Click the **Database Backup** button.



The Save As window appears.

	Save As	Search Backup	
	Organize ▼ New folder	¥≣ ▼ 🔞	—The File
	Recent Places	d Type Size 35 PM SQL Server Databa 7,616 KB	Directory drop- down menu
	Contenties Concents Contents		
The File— name field	Cecal Diak (C) SlowPoke (D) File name Seve as type: Database (*.mdf) Hide Folders	v Save Cancel	–The Save button

The Save As Window

- 4. Click the **File Directory** drop-down menu to navigate to the desired location where the backup file will be stored.
- 5. Enter the desired file name for the database backup file in the **File name** field.
- 6. Click Save.

The **Database backed up successfully** message appears.



The Database Backed Up Successfully Message

7. Click **OK** to return to the Database Backup/Restore screen.

Restoring the Database

Administrators can restore the simulator database using the Database Restore feature located on the Configuration screen.

To back up the simulator database:


1. From the icon dashboard, click the **Configuration** icon.



The Configuration Icon

2. From the Configuration screen, select the **Database Backup/Restore** option from the **Configuration** panel.

The Database Backup/Restore screen appears.



The Database Backup/Restore Screen

3. Click the **Restore Database** button.




The Open window appears.

GAE Healthcare - CathLa	abVR™ ▼					2/21/20	20 10:4	17 AM	⑦ Default Administrator • en •
*			1	D			×		
Con to per									
$e \rightarrow e \uparrow h$	is PC > Lecal Disk (C) >			νő	Search Local		2		
Organize - New folds	er.					81 • E	0		
Screenshots ^	Nama	Date modified	Type	Sax					
E Videos	CathLabVR - Edwards Lifesciences - Sapien		File folder						
c > CneDrive	latel	\$/23/2019 7/05 AM	File folder						
v 💷 This PC	Meditoric HIS PerfLops	11/14/2019 7/45 AM 3/18/2019 9/52 PM	File folder File folder						
> 30 Objects	Program Files	\$/23/2019 3:32 PM	File folder						
> Dektop	Program Files (c86)	12/11/2219 6:49 AM	File folder						
> R Decuments	Users	\$(20/2019 1:54 PM	File folder						
> 🕹 Dewnloads	Windows	12/13/2019 7.26 AM	Filefolder						
> Masic									
> E Pictures									
> Widees									
> 1_ Local Disk (C)									
> 💣 Netvolk									
> 🕐 Nebvork								Backup Database	
Fden	ane			Ŷ	Ostabase ("a Open	mdt) Can	cel		
								Restore Database	
	10000 C								

The Open File Window

- 4. Click the **File Directory** drop-down menu to navigate to the desired location where the backup file is stored.
- 5. Enter the desired file name for the restoration of the database in the **File name** field.
- 6. Click **Open**.

The **Back Up Current Data Warning** message appears.




The Back Up Current Data Warning Message

Note: The message states any data that has been added, deleted or modified since the last backup will be lost when a restore is performed.

Click **Cancel** to cancel the restore and return to the Database Backup/Restore screen to perform a backup.

7. Click **OK** to continue with the restore.

The **Database restored successfully** message appears.



The Database Restored Successfully Message

8. Click **OK** to return to the Database Backup/Restore screen.



Adjusting the Volume Controls

Administrators can configure the sound and volume hotkeys from the Configuration screen.

To adjust the volume controls:

1. From the icon dashboard, click the **Configuration** icon.



The Configuration Icon

2. From the Configuration screen, select the **Sound Volume Control** option in the **Configuration** panel.

2/23/2020 8:32:29 AM ⑦ Default Administrator • -0-Ċ Ð The Volume The Sound Increase Volume instructions Control option The Volume Decrease instructions The **Sound** Volume Hotkeys Configuration button

The Sound Volume Control screen appears.

The Sound Volume Control Screen

To increase the volume, hold down on the CTRL key and tap the up arrow key.

To decrease the volume, hold down on the CTRL key and tap the down arrow key.

To configure hotkeys for the sound volume:

1. Click the **Sound Volume Hotkeys Configuration** button.

The Sound Volume Hotkeys window appears.

	💜 Sound Volume Hotkeys		×	
The Enable Sound — Volume Indicator	✓ Enable Sound Volume Indicator	e indicator:		
checkbox	Effect during showing volume indicator:	Scroll effect	•	
	Effect during <u>h</u> iding volume indicator:	Scroll effect	•	
The Volume indicator — transparency radio button	Volume indicator transparency:	Previ	iew .	—The Preview button
	Vert Place Icon in the Taskbar			
The OK —	Active Keys - hotkeys software - dick	here for more info	<u>.</u>	
button	<u>A</u> bout OK	Cancel Hel	þ	

The Sound Volume Hotkeys Window

- 2. Select the **Enable Sound Volume Indicator** (selected by default) to display the volume increase/decrease slider when the hotkeys are used.
- 3. Select the Volume indicator transparency radio button.
- 4. Use the slider to adjust the transparency of the volume indicator when it appears on the screen.
- 5. Click the **Preview** button to see how the volume indicator will display with the adjusted settings.
- 6. Click **OK** to return to the Sound Volume Control screen.



Print

From the Print screen, administrators can print customized data including user profiles, usage summaries, and reports.

To use the printing functions in the software, a network connection must be established.


Print Preview

To create a print preview:

1. From the icon dashboard, click the **Print** icon.



The Print screen appears.



The Print Screen

2. Select the desired data to include in the printout by clicking in the associated checkbox.

The **Groups** options appear.

	CAE Healthcare - CathLabVR™ ▼	1/8/2014 12:39:52 PM	⑦ Default Administrator • en •
	🙏 🏚 🌣 📥	13	
The Data	Users & Groups Curriculum Configuration Print	Import/Export	
	Print Select the data to include in the printout		
options	User Profiles Summaries of Usage User Reports		
	Select the groups to inlcude		
The Groups_	All Groups		
options	Unassigned Users First Year Residents		
•	in risk real residents		
	PRINT PREVIEW:		
			Print

The Print Screen

3. From the **Groups** options, select the desired groups to include in the print preview or select **All Groups**.

The **Users** options appear.

4. From the **Users** options, select the desired users to include in the print preview or select **All Users**.

CAE Healthcare - CathLabVR™ ▼	1/8/2014 12:42:22 PM	⑦ Default Administrator * en *	
Liters & Groups Curriculum Configuration	Import/Diport		
Print			
Select the data to include in the printout			
Vser Profiles Summaries of Usage User Reports			
Select the groups to inlcude	Select the users to include	_	
I All Groups	All Users		
✓ Unassigned Users	III admin		——The Users
🗵 First Year Residents	guest JohnDoe		options
			options
PRINT PREVIEW:			
		<u>.</u>	
		Print	

The Print Screen



Note: If the User Reports data was selected, the **User Reports** options appear. From the **User Reports** options, select the desired user reports to include in the print preview or select **All User Reports**.

And Sources Consolution Consolution Period			
Oselect the data to include in the printout User Profiles □Summaries of Usage ℤUser Reports Oselect the groups to include ①Select the			
User Profiles □Summaries of Usage ℤUser Reports ②Select the groups to inkcude ③Select the			
③Select the groups to inicude ④Select the			
MARC	isers to include	Select the reports to include	
All Users		Z All Reports	The
		☑ JohnDoe: Sample Task - Case 1: 12/20/2013 6:28:13 AM	Daman
		▼ JohnDoe: Sample Task - Case 1: 12/19/2013 3:08:38 PM	Repor optior
Ø JohnDoe		☑ JohnDoe: Sample Task - Case 1: 12/19/2013 8:57:47 AM	ontion
State of the second sec		V JohnDoe: Sample Task - Case 1: 12/19/2013 8:52:56 AM	οριιοι
		JohnDoe: Sample Task - Case 8: 12/18/2013 4:24:28 PM	
		JohnDoe: Sample Task - Case 7: 12/18/2013 4:23:24 PM	
		JohnDoe: Sample Task - Case 6: 12/18/2013 4:22:36 PM	
		IohnDoe: Sample Task - Case 5: 12/18/2013 4:21:28 PM	
		JohnDon-Samola Tark Care A: 17/19/2012 A:20-24 DM	
RINT PREVIEW:			
ample Course 1 - Case 1: Mr. N - 61-year-old male :			
ohnDoe			
2/20/2013 6:28 AM			
Metrics			
SUMMARY Your Result			
Total procedure time: 0:00:43			
Complication occurred during procedure.: No			
ANGIOGRAPHY Your Result			
Aortography performed: No			1
CA angiography performed: No RCA angiography performed: No			1

The Print Screen

The selected data appears in the **Print Preview**.

CAE Healthcare - CathLabVR™ ▼		1/8/2014 12:45:05 PM	⑦ Default Administrator • en •
Contraction Contraction of	Configuration	to Import/Sport	
Print Select the data to include in the printout			
-			
Select the groups to inlcude		Select the users to include	Select the reports to include
Z All Groups		2 All Users	R All Reports
I Unassigned Users		🗷 admin	I JohnDoe: Sample Task - Case 1: 12/20/2013 6:28:13 AM
R First Year Residents		🛙 guest	Ø JohnDoe: Sample Task - Case 1: 12/19/2013 3:08:38 PM
		🗷 JohnDoe	JohnDoe: Sample Task - Case 1: 12/19/2013 8:57:47 AM
			✓ JohnDoe: Sample Task - Case 1: 12/19/2013 8:52:56 AM
			JohnDoe: Sample Task - Case 8: 12/18/2013 4:24:28 PM
			Ø JohnDoe: Sample Task - Case 7: 12/18/2013 4:23:24 PM
			Id JohnDoe: Sample Task - Case 6: 12/18/2013 4:22:36 PM
			JohnDoe: Sample Task - Case 5: 12/18/2013 4:21:28 PM
			JohnDoo: Samala Tark Case 4: 12/19/2012 4:20:24 Dk4
PRINT PREVIEW:			
Sample Course 1 - Case 1: Mr. N - 61-year	-old male :		
JohnDoe			
12/20/2013 6:28 AM			
Metrics			
metrics		_	
SUMMARY	Your Result		
Total procedure time: Complication occurred during procedure.:	0:00:43 No		
ANGIOGRAPHY	Your Result		
Aortography performed:	No		
LCA angiography performed:	No		
RCA angiography performed:	No		
			Print

The Print Screen

The document is ready to print.

Note: A network connection is required for printing documents from the simulator to a printer. To add a printer to the network, exit the software and consult your Windows Help menu.

Printing Data to a PDF File

If network connection is not available or no printers are available on the network, administrators can print data to a PDF file and save the file to an external device to print on a computer with a network printer connection.

To print data to a PDF file:

1. From the Print screen, click **Print**.

CAE Healthcare - CathLabVR™ ▼	1/8/2014 12:45:05 PM	⑦ Default Administrator • ended
Users & droups Curriculum Configuration	Pert Report/sport	
Print		
Select the data to include in the printout		
User Profiles 🔲 Summaries of Usage 📝 User Reports		
Select the groups to inlcude	Select the users to include	Select the reports to include
All Groups	Z All Users	I All Reports
I Unassigned Users	i admin	I JohnDoe: Sample Task - Case 1: 12/20/2013 6:28:13 AM
First Year Residents	₩ quest	JohnDoe: Sample Task - Case 1: 12/19/2013 3:08:38 PM
	⊯ JohnDoe	International State of the S
		☑ JohnDoe: Sample Task - Case 1: 12/19/2013 8:52:56 AM
		Id JohnDoe: Sample Task - Case 8: 12/18/2013 4:24:28 PM
		III JohnDoe: Sample Task - Case 7: 12/18/2013 4:23:24 PM
		Izi JohnDoe: Sample Task - Case 6: 12/18/2013 4:22:36 PM
		JohnDoe: Sample Task - Case 5: 12/18/2013 4:21:28 PM
		Jal. John Doo: Samela Tark Core A: 12/19/2012 A:20:24 Dk4
INT PREVIEW:		
mple Course 1 - Case 1: Mr. N - 61-year-old male :		
ohnDoe		
2/20/2013 6:28 AM		
Aetrics		
SUMMARY Your Result		
otal procedure time: 0:00:43		
omplication occurred during procedure.: No		
ANGIOGRAPHY Your Result Aortography performed: No		
CA angiography performed: No		

The Print Screen

The Print window appears.

9 Print					
General	Options				
Selec	t Printer				
	Fax				
0	Microsoft P	rint to PDF			
Statu Local Comm	tion:	leady		Print to file	Preferences Find Printer
	Range			Tarras - a	
• A		00.00	D	Number of copie	s: 1 ≑
	election	Current	Page		
	ages:	1		Collate	1 22 33
		le page number example, 5-12	or a single	1	ja <u>7</u> a 3a
				14	
			Pri	nt Cance	Apply

The Print Window



- 2. Select the Microsoft Print to PDF.
- 3. Click on the **Print** button.

The Save Print Output As window appears.

ZAE Save Print Output	As						×
← → ~ ↑ 💻	> This PC >			ٽ ~	Search This PC		,0
Organize 👻						₩. ₩.	0
 Medtronic HI: Music Screenshots Videos OneDrive This PC 30 Objects 	S Polders (7	7)) Objects wwnloads deos	Desktop		Documents Pictures		
 Job Objects Desktop Documents Downloads Music Pictures Videos 	Lo	nd drives (1) cal Disk (C:) 1 GB free of 232 GB					
File name:							~
Save as type:	PDF Document (*.pdf)					~
∧ Hide Folders					Save	Cancel	

The Save Print Output As Window

4. Select the desired folder where the file will be saved.

Note: If the location is on an external device, use the USB port located on the back of the simulator to connect the external device to the simulator.

5. Enter the desired file name in the **File name** field and click **Save**.

To open and view the saved PDF file:

- 1. After exiting the software, File Explorer may be used to navigate and locate the saved PDF file.
- 2. Double-click on the file name. If prompted to select the application to open the file with, select Adobe Acrobat.

The **Adobe Reader** window appears with the completed PDF.

Adobe Reader v Document Tools W	indow Help	-				
	/ 51 🕑 💌	1004		Find		
🧆 10 📥 🔸	/* 00	10075		rina		
						Page 1 of 51
	tions - Case 1: Mr.					
test :	bons - Case 1: NW.	N - 01-year-	ord mare :			
12/20/2013 6:28 AM	а					
Metrics						
SUMMARY			Result			
Total procedure tin		0:00:4	3			
ANGIOGRAPHY	red during procedu		Result			
Aortography perfor	and the second	Your	Result			
LCA angiography p	erformed:	No				
RCA angiography p	erformed:	No				
FLUOROSCOPY			Result			
Total fluoroscopy to Total cine time:	THE SHIT	0.00.0				
Total radiation exp	No. of Mr.	2.2 0				
Total amount of co	ntrast used:	0 mL				
Total # of stills save	d:	0				
CONTRAST LOG						
Contract was not used	during the intervention					
TOOL SELECTION	log					
No tools were selecte	d during the intervent	on.				
COMPLICATIONS	106					
No complications occ	arred during the inter-	ention.				
STENOSIS LOG						
Location	Dismeter Lengt	Bre Brit	Post-PCI	Lesion Coverage		
Proximal-Mid LAD		1 36.8 %	36.8 %	0.0 %		
Distal LAD	1.7 mm 15.2 m	13 32.9 %	32.9 %	0.0 %		
Distal LCx	2.0 mm 20.3 m		32.5 %	0.0 %		
Proximal RCA Mid RCA	3.9 mm 11.0 m	m 39.0 %	39.0 %	0.0 %		
ADJUNCT THERAS					-	
	cured during the inte	unition.				
IMAGES	and an and an and					
	red during the proces					
PROCEDURE LOG	the party are proces					
Advent Recording And						
00:07 CINE: started 00:10 CINE: Stopped	RR 5. 0802 355. At					
00:15 VITALS: NR 47 00:01 CIME: started	RR 9, 0902 999, M	110/01(94) #	with g			

The Adobe Reader Window

Printing Data to a Network Printer

To print data to a network printer:

1. From the Print screen, click the **Print** button.

CAE Healthcare - CathLabVR™ ▼	1/8/2014 12:45:05 PM	⑦ Default Administrator * en
Lisers & Groups Curriculum Configuration	Pire Report Report	
Print		
Select the data to include in the printout		
🗉 User Profiles 🛛 Summaries of Usage 🗷 User Reports		
Select the groups to inlcude	Select the users to include	Select the reports to include
I All Groups	2 All Users	2 All Reports
軍 Unasigned Users ※ First Year Residents	iii admin iii guest iii JohnDoe	IdehnDoe: Sample Task - Case 1: 12/20/2013 6/28:13 AM IdehnDoe: Sample Task - Case 1: 12/29/2013 3/8:38 PM IdehnDoe: Sample Task - Case 1: 12/29/2013 3/8:374 AM IdehnDoe: Sample Task - Case 1: 12/29/2013 8/5256 AM IdehnDoe: Sample Task - Case 1: 12/29/2013 8/5256 AM IdehnDoe: Sample Task - Case 1: 12/29/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM IdehnDoe: Sample Task - Case 1: 12/28/2013 4/22/38 PM
PRINT PREVIEW:		
Sample Course 1 - Case 1: Mr. N - 61-year-old male :		
JohnDoe 12/20/2013 6:28 AM Metrics		
SUMMARY Your Result		
Total procedure time: 0:00:43		
Complication occurred during procedure.: No		
ANGIOGRAPHY Your Result Aortography performed: No.		
LCA angiography performed: No		

The Print Screen



The Print window appears.

eneral Options			
Select Printer			
😸 Printer 1			
PDF Complete			
•		III	Þ
Status: Ready		Print to file	Preferences
Location:			
Comment:			Find Printer
Page Range			
All		Number of copies:	1 🚔
Selection	Current Page		
O Pages:	1	Collate	3-63-63
Enter either a single pag page range. For examp		1	1 22 33

The Print Window

- 2. Select the desired network printer.
- 3. Click **Print**.

The document prints to the network printer.

Note: A network connection is required for printing documents from the simulator to a printer.

Import/Export

The Import/Export tab gives access to import/export capabilities for handling user data and simulation results.

Importing Users and Groups

Using the import function, administrators can add data for multiple users in one process.


To import users:

1. From the icon dashboard, click the **Import/Export** icon.



The Import/Export Icon

The Import/Export screen appears.



The Import/Export Users Screen

2. From the Import/Export screen, click Users & Groups.



The **Import** and **Export** buttons appear.

	CAE Healthcare - CathLabVR [™] ▼	1/8/2014 12:57:19 PM	② Default Administrator • en •
	Lien à Geogra Carriculum Configuration Print	Provent Street	
	Select the elements you want to import or export		
		USERS & GROUPS RESULTS (Export)	
	Select if you want to import or export users & groups		
The Import—	I	IMPORT EXPORT	
button			
batton			

The Import/Export Users Screen

3. Click Import.

The **File Selection** tools appear.

	CAE Healthcare - CathLabVR [™] ▼	1/8/2014 12:58:53 PM	⑦ Default Administrator • en •
	💵 📾 🌣	ا ال	
	Users & Groups Curriculum Configuration	First Import/Expert	
	Select the elements you want to import or export		
		USERS & GROUPS RESULTS (Export)	
	Select if you want to import or export users & groups		
		IMPORT EXPORT	
The Durance	Select a file to import users from		
The Browse—	Import users from:		
button	Create a new group from this list	wse	
	Import		

The File Selection Tools

4. Click **Browse**.


The **Select File** window appears.

5. Select the comma delineated (CSV) file to import.

Note: If the file is located on an external device, use the USB port located on the back of the simulator to connect the external device to the simulator.

6. Click **Open.**


The file name appears in the **Import users from** field.



The File Selection Tools

- 7. Select **Create a new group from this list** to create new group from the file, if desired.
- 8. Click Import.

The Import Success message appears.



The Import Success Message

9. Click **OK.**

The users are imported into the software.





Exporting Users and Groups

Administrators can export data to an external device for evaluation and record purposes. Data for individual users or groups can be customized and exported into a comma delineated (CSV) spreadsheet file from the CathLabVR simulator.

To export users and groups:

1. From the icon dashboard, click the **Import/Export** icon.



The Import/Export Icon

The Import/Export screen appears.

	CAE Healthcare - CathLabVR™ ▼			1/8/2014 12:51:45 PM	@ D	efault Administrator 🔻 🔄 🔻
	Users & Groups Curriculum	Contiguation Print	inguest/Expect			
The Users &—	Select the elements you want to imp	port or export				
Groups button				USERS & GROUPS RESULTS (Export)		

The Import/Export Screen

Using the Administrator Tools

2. From the Import/Export screen, click **Users & Groups**.

	CAE Healthcare - Ca	athLabVR™ ▼				1/8/2014 12:57:19 PM	1	⑦ Default Administrator	r 🔻 en 🔻	
	Users & Groups	Curriculum	Configuration	.	Import/Expert					
The Users &—	Select the element	nts you want to im	port or export			USERS & GROUPS RESULTS	Funget			
Groups button	Select if you w	vant to import or e	eport users & groups			USENS & GROUPS RESULTS	export			
						IMPORT EXPORT				–The Export
										button

The Import/Export Users Screen

3. Click Export.

The **Groups** options appear.

	CAE Healthcare - CathLabVR™ ▼	1/8/2014 1:04:02 PM	⑦ Default Administrator • en •
	🙏 👜 🌣 📥	1	
	Uvers & Groups Curriculum Configuration Print Select the elements you want to import or export	Import/Expert	
		USERS & GROUPS RESULTS (Export)	
	③Select if you want to import or export users & groups		
The Groups options		IMPORT EXPORT	
options	Select the groups to export		
·	Dunassigned Users		
	First Year Residents		

The Export User Screen

4. Select the group or **All Groups** from the **Groups** options.



The **Users** options appear.

	CAE Healthcare - CathLabVR™ ▼	1/8/2014 1:05:00 PM	② Default Administrator • en •
	Uura la Despa	Impositioned	
	Select the elements you want to import or export		
		USERS & GROUPS RESULTS (Export)	
	③Select if you want to import or export users & groups		
The Users —		IMPORT EXPORT	
options	③Select the groups to export	Select the users to export	
options	All Groups June 2019 June 201	All Users JohnDoe	
	⊯ First Year Residents	🗈 admin	
		🗉 guest	

The Export User Screen

5. Select the user or **All Users** from the users options.

The **File Selection** tools appear.

	CAE Healthcare - CathLabVR [™] ▼	1/8/2014 1:05:26 PM	⑦ Default Administrator • en •
	Uner & Girages Carriculum Carriguesion Print	Depart Growt	
	Select the elements you want to import or export		
		USERS & GROUPS RESULTS (Export)	
	Select if you want to import or export users & groups		
		IMPORT EXPORT	
	③Select the groups to export	Select the users to export	
	I All Groups		
	Unassigned Users	I JohnDoe	
	First Year Residents	Ø admin	
		Ø guest	
	Select a file to export users to		
The Browse —	Export users to:		
	Browse		
button	Export		

The User Export Screen

6. Click **Browse**.

The **Select File** window appears.

- 7. Type the name of the file in the **File name** field.
- 8. Select the location where the file will be saved when the exporting is complete.

Note: If the location is on an external device, use the USB port located on the back of the simulator to connect the external device to the simulator.

9. Click Save.


The file name appears in the **Export users to** field.

	CAE Healthcare - CathLabVR™ ▼	1/8/2014 1:05:26 PM	② Default Administrator • en •
	Jure à Geogra Carriculari. Configuration Print	Lingert/Tippert	
	Select the elements you want to import or export		
		USERS & GROUPS RESULTS (Export)	
	OSelect if you want to import or export users & groups	IMPORT EXPORT	
	③Select the groups to export	③Select the users to export	
	IF All Groups	2 All Users	
	Unassigned Users First Year Residents	I JohnDoe I admin	
	R First Year Kesidents	i≪ admin ⊮ guest	
	Select a file to export users to		
	Export users to:		
The Export — button	Brows		

The Selected Export File

10. Click Export.

The **Export Success** message appears.



The Export Success Message

11. Click **OK**.

The file is exported.



Exporting Results

Administrators can export results to an external device for evaluation and record purposes. The data is exported into a comma delineated (CSV) spreadsheet file from the CathLabVR simulator.

To export results:

1. From the icon dashboard, click the **Import/Export** icon.



The Import/Export Icon

CAE Healthcare - CathLabVR [™] ▼	1/8/2014 12:51:45 PM	⑦ Default Administrator • en •	
11 💼 🌣 📥	1		
Contraction Contraction Contraction Pro-	USERS & GROUPS RESULTS (Suport)		—The Results (Export) button

The Import/Export Screen

Using the Administrator Tools

2. From the Import/Export screen, click **Results (Export).**

	CAE Healthcare - CathLabVR [™] ▼	1/8/2014 1:10.52 PM	⑦ Default Administrator • en •	
The Export Type — drop-down menu	Select the elements you want to import or export Oselect the users lat's export type Operating the elements of	USERS & GROUPS RESULTS (Supur)		—The Results (Export) button
The Groups — options				

The Export Results Screen

- 3. Select the export type from the drop-down menu. Choose from **Clear** (default), **Anonymous** or **Double Blind**.
- 4. Select the group or **All Groups** from the **Groups** options.
- 5. Select a user or **All Users** from the **Users** options.


	CAE Healthcare - CathLabVR [™] •	1/8/2014 1:13:07 PM	⑦ Default Administrator • en •
	Umr A Groups Curriculum Configuration Print	Line and the set	
	Select the elements you want to import or export		
		USERS & GROUPS RESULTS (Export)	
	Select the users list's export type		
The Users –	Export Type: Clear		
	Select the groups to export	Select the users to export	
options	If All Groups	2 All Users	
	Unassigned Users	I InnDoe	
	III First Year Residents	i i admin	
		If guest	
The Courses —	- + -1		
options	Select the courses to export		
options	ID All Courses		
	Sample Course 1		
	Sample Course 2		
	Unassigned Tasks		
	Sample Course 3		
	Sample Course 4		
	Sample Course 5		

The Export Results Screen

6. Select a course or **All Courses** from the **Courses** options.



7. Select a task or **All Tasks** from the **Tasks** options.



The Export Results Screen

- 8. Click Browse.
- 9. Select the location where the file will be saved when the exporting is complete.

Note: If the location is on an external device, use the USB port located on the back of the simulator to connect the external device to the simulator.

10. Click Open.


The file name appears in the **Export results to** field.

Healthcare - CathLabVR [™] ▼	1/8/2014 1:13:33 PM	② Default Administrator • en •	
Select the elements you want to import or export	USERS & GROUPS RESULTS (Export)		
Select the users list's export type			
Export Type: Clear			
Select the groups to export	Select the users to export		
R All Groups	R All Users		
Unassigned Users	I JohnDoe		
First Year Residents	🗷 admin		
	2 guest		
		Select a folder to export results to	
		Export results to:	
		C:\Users\cae\Documents Browse	
		The Ex	nor
		ExportThe Ex	hou
Select the courses to export	Select the tasks to export	button	1
R All Courses	□ All Tasks	Dutton	
	😿 Sample Course 1 Task 1 - Case 1	*	
☑ Sample Course 2	Sample Course 1 Task 1 - Case 2	3	
Unassigned Tasks	☑ Sample Course 1 Task 1 - Case 3		
☑ Sample Course 3	Izi Sample Course 1 Task 1 - Case 4		
☑ Sample Course 4	🗷 Sample Course 1 Task 1 - Case 5		
I Sample Course 5	🗵 Sample Course 1 Task 2 - Case 1		
	📧 Sample Course 1 Task 2 - Case 2		
	😿 Sample Course 1 Task 2 - Case 3		
	- County County 1 Tech 3, County	*	

The Export Results Screen

Using the Administrator Tools

11. Click Export.



The Result Export Success Message

12. Click **OK.**



Exiting the Software

Administrators should always log out of the software once they are finished using the software.

To exit the software:

1. From any screen, click on the username in the top right corner of the screen.

The **Logout** drop-down menu appears.



The Home Screen with Logout Drop-Down Menu

2. Select **Logout** from the **Logout** drop-down menu.

The administrator is logged out of the software and the Login screen appears.

en CAE CathLab \R series
Username
Password
Login
Copyright @ CAE Healthcare Inc., 2019. All Rights Reserved. 1.0.13.0

The Login Screen

THIS PAGE INTENTIONALLY LEFT BLANK





APPENDIX A - FAQ

This Frequently Asked Questions (FAQ) section addresses some common concerns and inquiries. For any questions not addressed here, please contact CAE Healthcare Customer Service. Contact information is provided on the last page of this document.

Dual Monitor Configuration

When a secondary monitor is connected to the laptop, the software defaults to dual monitor configuration. If the dual monitor software is not set up, both monitors have the same image.

To configure the monitors:

- 1. Click the Windows Notification button on the lower-right corner of the monitor.
- 2. Click on Project.
- 3. If you do not see the Project option, click **Expand** to view more options.

Ę	ä	ŵ	Ľ
Tablet mode	Network	All settings	
<i>ф</i>		৶	
Airplane mode		Alarms only	
*	¢.	*	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~
50%	Night light	Not connected	VPN
φ	Ş	Ð	г ө
Battery saver	Screen snip	Project	Connect
名 ⁹ へ 幅 ╔ ロゥ) ^{4:51 PM} 🗔			

Windows Notifications and Settings

a. Click on **Extend**.





For more information about CAE products, contact your regional sales manager or the CAE distributor in your country, or visit caehealthcare.com. Tel +1 941-377-5562 or 866-233-6384

For customer service, please contact CAE.

Customer Service Headquarters - United States

Monday - Friday from 7:00 a.m. to 6:00 p.m. ET Phone 1-866-462-7920 Email: srgcustomerservice@cae.com

Customer Service - Canada Monday - Friday from 8:00 a.m. to 5:00 p.m. ET Phone 1-877-223-6273 Email: can.service@cae.com

Customer Service - Europe, Middle East, and Africa

Monday - Friday from 8:00 a.m. to 5:00 p.m. CET Phone +49 (0) 6131 4950354 Email: international.service@cae.com

Customer Service - United Kingdom and Ireland

Monday - Friday from 9:00 a.m. to 5:00 p.m. GMT Phone +44 (0)800-917-1851 Email: uk.service@cae.com

Customer Service - Latin America

Monday - Friday from 9:00 a.m. to 5:00 p.m. BRT/BRST Phone +55 11 5069-1510 Email: la.service@cae.com

Customer Service - Asia Pacific

Monday - Friday from 8:00 a.m. to 5:00 p.m. CET Phone +49 (0) 6131 4950354 Email: ap.service@cae.com

©2020 CAE

