

## Recommended sizes for clinical supplies

Item	CAE Apollo Nursing	CAE Apollo Prehospital	CAE Ares	CAE Aria	CAE Blue Phantom	CAE HPS	CAE Juno	CAE Lucina	CAE Luna	CAE Medicor
Urinary Catheter	16 Fr	16 Fr	16 Fr	10 Fr		14-16 Fr	16 Fr	14-16 Fr	6.5 Fr	14 Fr
Nasogastric Tube	14 Fr*	14 Fr**	14 Fr**	10 Fr**		14 Fr**	14 Fr*	14 Fr**	8 Fr	
Gastric Lavage/ Gavage	14 Fr*						14 Fr*			
Airway Suctioning	14 Fr*						14 Fr*			
ETT	7.0-7.5 mm	7.0-7.5 mm	6.5-7.5 mm	5.5-6.0 mm		7.0-8.0 mm	7.0 mm	6.5-7.5 mm	3.5 mm	
LMA	#4	#4	#4, #5	#2.5		#3		#4, #5	#1	
King Airway	#4	#4	#4	#3				#4		
Combitube®	37 Fr	37 Fr	41 Fr			37 Fr		41 Fr		
i-gel® Supraglottic Airway	#4	#4	#4	#2.5		#3		#4		
Oropharyngeal Airway	90 mm	90 mm	100 mm	70 mm		90 mm	100 mm	90 mm	50 mm	
Nasopharyngeal Airway	30 Fr	30 Fr	28 Fr	24 Fr		30 Fr	28 Fr	28 Fr	14 Fr	
Tracheostomy Tube	6 mm	6 mm	6 mm	5 mm		6 mm	6 mm		3.5 mm	
IV Cannula	18-22 gauge	18-22 gauge	20-22 gauge	22 gauge		18-22 gauge	20-22 gauge	18-22 gauge	24 gauge	22-gauge IV needles
IO Needle Tibial IO Needle Humerus		25 mm 25 mm	N/A 25 mm	15 mm 15 mm						
Chest Tube	28 Fr	28 Fr		20 Fr		24-28 Fr	28 Fr		8 Fr	
Needle Decompression	14 gauge 6 cm	14 gauge 6 cm	14 gauge 6 cm	16 gauge 3 cm	18-21 gauge	14 gauge 6 cm			22-24 gauge IV catheter	
Pericardiocentesis						16-18 gauge 10-15 cm				
Intrauterine Balloon								24 Fr		
Umbilical Catheter									5 Fr	
Gastrostomy Feeding Tube (PEG)				MIC-KEY* 0124-12						

\* With Fluid Return \*\* Insertion Only