

Recommended Sizes for Clinical Supplies

Item	CAEiStan™	CAEPediaSIM™	CAEBabySIM™
Urinary Catheter	16 Fr	10 Fr	
Nasogastric Tube Insertion	14 Fr **	10 Fr **	8 Fr **
Gastric Lavage / Gavage			
Airway Suctioning			
ETT	7.0 - 7.5 mm	5 mm Uncuffed	3.5 mm Uncuffed
LMA Unique	#4	#2	#1.5
King LTS-D LT-D	#4		
Combitube	37 Fr		
i-gel Supraglottic Airway	#4		
Oropharyngeal Airway	90 mm	60 mm	43 mm
Nasal-Pharyngeal Airway	30 Fr 7.5 mm	24 Fr 6 mm	20 Fr 5 mm
Tracheostomy Tube	6 mm	3.5 mm	
IV Cannula		20 - 22 g	20 - 22 g
Chest Tube	28 Fr	20 Fr	12 Fr
Needle Decompression	14 g 6 cm	14 g 6 cm	14 g 6 cm
Pericardiocentesis			
Intrauterine Balloon			
Umbilical Catheter			

* With fluid return

**Insertion only