

## Recommended Fluid Amounts for Reservoirs/Tanks

Reservoir	iStan®	Apollo® PreHospital	Apollo® Nursing	Caesar™	Lucina™	Athena™	Juno™
Airway			40 mL				30 mL
Gastric			60 mL				150 mL
Urinary					70 mL	70 mL	100 mL fluid 80 mL air
Ostomy							200 mL
Trauma Reservoir				1.4 L			
Right Leg Trauma	1.8 L	1.5 L					
Left leg Trauma	0.8 L						
IV Prime	60 mL	60 mL	60 mL	60 mL	50 mL	50 mL	
Hemorrhage Tank					1.8 L		