



CAE Blue Phantom

General Pathology Transvaginal Ultrasound Training Model

Realistic and resilient ultrasound training models

CAE Healthcare's Blue Phantom ultrasound training phantoms are built with realism and durability so you can count on their lifespan. The Blue Phantom trainers support at least 20 medical specialties including central line insertion, breast biopsy, and transvaginal ultrasound exams. CAE's realistic ultrasound training phantoms realistically replicate human tissue in every way. Built with a patented self-sealing Simulex™ tissue, the training models image like real tissue using any clinical ultrasound system and cannulate like real tissue to support hands-on training. All training models, except Endovaginal and FAST, are available in a transparent configuration, upon customer request.

For more information on Blue Phantom ultrasound training models visit bluephantom.com.

Your worldwide training partner of choice



What's in the Box

- General Pathology Transvaginal Ultrasound Model
- Soft Storage Case with Foam Block
- Privacy Skirt



General Pathology Transvaginal Ultrasound Training Model



- Extremely realistic imaging characteristics
- Constructed using Blue Phantom's patented ultra-realistic simulated human tissue providing users an extremely realistic training model
- Feels like real human tissue
- Anatomically accurate external anatomy
- Endovaginal and cervical canal
- Bladder with bladder wall
- Uterus with prominent endometrium
- Multiple uterine fibroids of varying sizes and sonographic characteristics
- Free fluid in cul-de-sac
- Left and right ovaries
- Ovarian follicles
- Ovarian cyst
- Large ovarian mass/abscess
- Broad ligaments
- Bowel and colon
- Accessory structures
- Extremely durable; use for repeated training
- Easy to clean – though not routinely necessary
- No replacement parts necessary
- Soft storage case and privacy skirt included
- Size 18" x 18" x 12" (46cm x 46cm x 30cm) (L x W x H)
- Weight 22 lbs. (10 Kg)